



The adjustable resin line clamp enables precise control or complete stoppage of individual resin flow lines during infusion processes. Compatible with polyethylene, nylon, and silicone tubes, it allows exact resin dosing.



TECHNICAL DATA	VALUE
Material	Stainless steel / Galvanized steel / POM
Flange	20 mm
Max. diameter of tube	16 mm

**STORAGE**

Shelf life: unlimited.

Storage conditions: store at 10–30°C in original packaging

**NOTES**

Technical values are provided to the best of our knowledge and are based on data considered reliable. Users are responsible for verifying suitability and assume all associated risks.